




The optimal choice for ANA diagnostics

Our wide range of screening and confirmatory tests



Screening test on gold standard level

- **Indirect immunofluorescence test (IIFT)** with human epithelial cells (HEp-2) for the differentiation of 15 nuclear, 9 cytoplasmic and 5 mitotic patterns according to ICAP *
- Simplified evaluation of IIFT patterns on HEp-2 cells due to their combination with tissue sections in BIOCHIP mosaics
- More mitotic cells due to the use of HEp-20-10 cells (HEp-2 cell variant with an increased division rate)

* International consensus on ANA pattern: www.anapatterns.org



Monospecific confirmatory tests from the full range of methods

ELISA & ChLIA

- Antigens in a mix or in separate profiles for the efficient detection of the most prevalent antibodies
- Single-parameter tests for quantitative determination

EUROLINE

- Multiplex tests for the parallel determination of up to 23 autoantibodies in one incubation
- Large choice of parameter combinations: from comprehensive ENA profiles to profiles for specific indications, such as for myositis or systemic sclerosis diagnostics



Flexible automation and software solutions for each test system

For more information and products, visit
<https://www.euroimmun.com/products/autoimmune-diagnostics/>

